

Cheetah Top Load Balances

YP Series

General Purpose Laboratory Balance

- Precision: 10 or 100 mg
- Range: 0—6,000 grams (refer to table below)
- Pan Sizes: refer to table below
- Overload protection
- Auto external calibration

CHEETAH

Reliability, Precision and Usability

Model	YP102N	YP202N	YP402N	YP601N	YP1201N	YP502N	YP802N	YP1002N	YP2001N	YP3001NJ	YP6001N
Precision (mg)	10			100		10			100		
Range (g)	0-100	0-200	0-400	0-800	0-1000	0-500	0-800	0-1000	0-2000	0-3000	0-6000
Pan Size (mm)	Ø110			Ø130		130 x 130			160 x 170		
Price + GST	\$332 + \$33.20	\$395 + \$39.50	\$457 + \$45.70	\$240 + \$24.00	\$271 + \$27.10	\$364 + \$36.40	\$395 + \$39.50	\$457 + \$45.70	\$302 + \$30.20	\$332 + \$33.20	\$395 + \$39.50

Cheetah Analytical Balances

FA Series

Great value for 4 decimal place accuracy

- Precision: 0.1 mg
- Range: 0-120 or 0-220 grams
- Pan size: Ø80 mm
- Large LCD display
- One touch operation
- Overload Protection
- Auto calibration (external weight)
- Level Indicator

Model	FA1204B	FA2204B
Precision (mg)	0.1 mg	
Range (g)	0-120	0-220
Pan Size (mm)	Ø80	
Repeatability (mg)	±0.1	
Linearity (mg)	±0.2	
Price + GST	\$1529 + \$152.90	\$1645 + \$164.50

JA Series

When 3 decimal place accuracy is sufficient

- Precision: 1 mg
- Range: 0-260 or 0-500 grams
- Pan Size: Ø80 mm
- Large LCD Display
- One touch operation
- Overload protection
- Auto calibration (external weight)
- Level Indicator

Model	JA2603B	JA5003B
Precision (mg)	1 mg	
Range (g)	0-260	0-500
Pan Size (mm)	Ø80 mm	
Repeatability (mg)	±1	
Linearity (mg)	±2	
Price + GST	\$1111 + \$111.10	\$1520 + \$152.00

Precisa Top Load Balance

165 BJ Series

Supports Animal Weighing

- Precision: 1 mg to 1 g
- Range: 0-12100 g
- Pan size: refer to table below
- Backlit LCD display
- Supports IQ/OQ/PG
- Animal weighing
- Swiss technology
- Loads of other features



Model/ Code	Capacity	Readability	Repeatability	Linearity	Pan size (mm/inch)	Price + GST
BJ100M	102 g	1 mg	1 mg	2 mg	Ø100 / 3.9	\$750 + \$75
BJ410C	410 g	0.01 g	0.01 g	0.02 g	135 x 135	\$620 + \$62
BJ010C	610 g	0.01 g	0.01 g	0.02 g	135 x 135	\$660 + \$66
BJ1000C	1020 g	0.01 g	0.01 g	0.02 g	135 x 135	\$710 + \$71
BJ1200C	1200 g	0.01 g	0.01 g	0.02 g	135 x 135	\$750 + \$75
BJ2200C	2200 g	0.01 g	0.01 g	0.02 g	135 x 135	\$800 + \$80
BJ2100D	2100 g	0.1 g	0.1 g	0.2 g	170 x 140	\$590 + \$59
BJ4100D	4100 g	0.1 g	0.1 g	0.2 g	170 x 140	\$620 + \$62
BJ6100D	6100 g	0.1 g	0.1 g	0.2 g	170 x 140	\$660 + \$66
BJ8100D	8100 g	0.1 g	0.1 g	0.2 g	170 x 140	\$710 + \$71
BJ6100G	6100 g	1 g	1 g	2 g	170 x 140	\$590 + \$59
BJ12100G	12100 g	1 g	1 g	2 g	170 x 140	\$620 + \$62

Density application optional

Precisa Analytical Balance

360 ES Series

*D decimal Place Accuracy
Ideal for Pipette Calibration*

- Precision: 0.01 mg
- Range: 125 or 102/225 grams
- Pan Size: refer to table below
- Supports IQ/OQ/PQ
- Animal weighing
- Swiss Made
- Loads of features and accessories



Model	Capacity	Readability	Repeatability	Linearity	Pan Size mm/inch	Response time	Price + GST
ES-125SM	125 g	0.01 mg	0.03 mg	0.06 mg	Ø80 / 3.1	10 s	\$4000 + \$400
ES 225SM-DR	102 g / 225 g	0.01 mg / 0.1 mg	0.03 / 0.1 mg	0.06 / 0.2 mg	Ø80 / 3.1	10 s / 3 s	\$3600 + \$360